Quick Guide for TEM

Dr. Fritz Schulze-Wischeler

Content

1. Preliminary Remarks.....	2
2. Occupational Safety.....	2
TEM Bright Field.....	3
3. Starting the Session.....	3
4. Place the specimen in the single-tilt specimen holder	6
5. Probenhalter in die Compustage einbringen.....	7
6. Control Panel.....	9
7. Alignments.....	10
8. Messen mit Digital Micrograph.....	17
9. Contrast Aperture.....	19
10. Ending the Session.....	20
TEM other Modes.....	23
11. Dark-Field Imaging.....	23
12. Selected Area Diffraction	24
13. STEM.....	26
14. STEM-EDX.....	29
Troubleshooting & Maintenance.....	33
15. Known Problems & Solutions.....	33
16. Restart TEM Server.....	34
17. Restart FEG and High Tension.....	35
18. Set up new TEM user on TEM control computer.....	37
19. Recalibrate Camera References.....	41
20. Change of an Aperture.....	43
21. Maintenance.....	44
More Information.....	45
22. Contact information.....	45
23. Features of the TEM	45
24. Book a Session.....	45
25. Analysis Software	47

1. Preliminary Remarks

- Only instructed users may work on the TEM.
- This quick guide is intended as a quick reference for trained users. Press F1 to access the Online Help on the TEM computer.
- The instructions on occupational safety must be fully understood.
- Please do not eat or drink at the TEM (it is ok in the anteroom).
- If you want to leave the TEM room, the column valves must *always* be closed.
- The session must be booked and the cost coverage by the institute/working group must be clarified in advance. (220 Euro per session, plus sample preparation if necessary).
- The instructions for publications must be known and followed. Usage regulations: "When using the TEM, users agree to thank the LNQE in publications."
- Things that can be damaged most quickly on the TEM due to incorrect operation:
 - Bending the holder - Therefore, always handle the holder carefully and never with force.
 - Vacuum collapse on the FEG - Therefore, always keep Column Valves closed when changing anything on the vacuum system.
 - Red light on the Compustage means: Never insert or remove a specimen holder in the Compustage.
 - Damage scintillator on camera due to too much beam power - Therefore, double check that the shield is closed when switching to Diffraction Mode.

2. Occupational Safety


- **X-ray protection:** The initial briefing for LNQE equipment subject to radiation or X-ray protection is carried out by the technical manager. Thereafter, the instruction must be refreshed to the day at the latest. These briefings are then possible online.
- **Liquid nitrogen LN2:** briefing by Fritz Schulze-Wischeler. The hazards are: Cold burns, lack of oxygen due to displacement, explosion hazard for closed containers, O₂ accumulation on cold surfaces. => Wear protective cryogenic gloves and goggles when filling.
- **SF6:** There is SF6 in the TEM at the FEG and in the equipment room. SF6 has two hazards:
 - Risk of suffocation: If it escapes, it displaces oxygen and there is a risk of suffocation (SF6 is heavier than air).
 - Fluorine poisoning: SF6 decomposes to sulfur and fluorine at over 250 °C, releasing highly reactive fluorine. Fluorine is very toxic.
- **In case of emergency:** In case of fire or other emergency situation (building evacuation), press the OFF key on the System On/Off panel. The microscope will shut down, minimizing any danger to persons in its vicinity or to yourself. Switch off the TEM PC. Immediately leave the room with the microscope.


TEM Bright Field

3. Starting the Session


- **Check that the HT and Vac switches on the control panel are illuminated**
This only works if the room light is on. If not both control lamps are lit, do not continue, but contact the TEM instrument officer.

- **Log in to the computer in the anteroom and check the asset ledger**
Logging into the LNQE domain. In the "TEM" group drive, the asset ledger "Anlagenbuch TEM FEI Tecnai F20.XLS" is located in the main directory. Check here (especially the "Special" column) in which state the TEM is currently is.
- **Log in to the TEM computer and start control programs**
 1. log in with the individual user name.
 2. start programs (in this order):
 1. TEM User interface
 2. digital micrograph (wait a while until the menus are loaded)
 3. TIA
- **Check whether "Column Valves" are closed, HT and FEG are on**

- Transfer important values to the asset ledger.


Write the vacuum values in the asset ledger.

Write the FEG emission current in the asset ledger. If the current has changed significantly from the last session, contact the TEM instrument officer.

- Check apertures

C2 to 4 for TEM, Objective & Selected Area aperture **out**


Condenser Aperture 1 (C1, always at position 4)

Condenser Aperture 2 (C2 always on position 3).

Objective Aperture

Selected Area Aperture


Aperture is:

Inside

Out

- Check Compustage

The Compustage must be in the zero position, i.e. X, Y, Z and α , β must all be approximately 0. If this is not fulfilled, then a Holder reset must be performed (see under "Ending the session").

- **Fill liquid nitrogen**

1. Cover the viewing screen on the TEM with the rubber cover.
2. Put on protective goggles and gloves.
3. Carefully fill the cold trap on the TEM to approx. 1 cm below the edge.
4. After approx. 10 minutes, you can load the trap.

- **Open vacuum overview**

(At the bottom right of the screen)

The different values represent logarithmic values of pressures.

IGP = Ion getter pump

Pen = Penning vacuum gauge

Pir = Pirani vacuum gauge

IGP2 and IGP3 are the FEG vacuum (IGP2 must be 1)


IGP1 is the column vacuum (must be 6 if ideally pumped down)

V7 + V4 closed correspond to the "Column Valves: Closed" state

P3 is the vacuum at the camera


P1 is the vacuum at the buffer tank

P2 is the backing line vacuum

Open the vacuum overview here


4. Place the specimen in the single-tilt specimen holder


- The brass-colored part between the O-ring and the head end of the holder must never be touched with the hand or other parts that are not suitable for vacuum. wear protective gloves if necessary.
- Press lightly against the black end of the holder with one hand to prevent it from slipping out of the protective tube.
- Insert the tool (located in the holder, see illustration) into the hole at the head end of the clamp. Then **carefully and without force**, lift the clamp straight up until it is at right angles to the holder (see photo).
- Place the specimen as centered as possible in the circular depression of the tip of the specimen holder. Lift the sample with the small sample suction pump as close to the edge as possible (or alternatively with tweezers). For TEM & STEM, the sample should be placed around so that the material to be microscoped is now facing down. When mounting in the TEM, the sample holder is rotated 180° so that the material is now facing upwards.
- Lower the clamp with the tool **carefully and without force** straight down onto the sample. Guide the clamp during this process and do not allow it to snap. Make sure that the specimen remains correctly in position. The spring in the specimen holder presses the clamp against the specimen, so it is not necessary to press it yourself. Put the tool back into the holder.
- Gently pull the holder out of the protective tube and rotate the holder 180° so that the specimen is facing down. Tap the black end of the holder a few times. Turn the holder back and check that the specimen has not moved (movement is an indication that it was not properly inserted and clamped).
- Pull the holder out of the protective tube and inspect the O-ring closely for dirt, hair, etc. under the magnifying glass. Dirt can be carefully removed with a tweezer or toothpick.
- If the O-ring has become too dry: Never grease the O-ring independently, but contact the TEM instrument officer.

5. Insert specimen holder into the Compustage

- Touch the Selected Area control once with your left hand to discharge any static charges. Now stand still.
- Carefully insert the specimen holder into the compustage with the small pin near the holder tip at the 5 o'clock position so that it fits into the small notch inside the compustage.
- Alternatively: Carefully insert the specimen holder with pin at 6 o'clock position until you feel a resistance at the end. Then turn to 5 o'clock so that the specimen holder can be inserted all the way in. (Here the warning light is on briefly - this is ok and normal).
- The pump starts automatically and evacuates the Compustage. The red warning lamp is on: Do not move the specimen holder! The Single-Tilt-Holder can be released carefully. The double-tilt holder should be held carefully for the entire pumping time.
- The respective holder (single or double) must be selected on the control interface on the PC and confirmed with the small arrow button ↵.


- For Double-Tilt only: Connect the cable from the holder to the Compustage. The cable must not block later or touch the protective cap (otherwise vibrations will be coupled in).


- The remaining **pumping time** is displayed at the bottom of the vacuum override (**180 seconds**). (Experienced users set "save numbered" in Digital Micrograph during this time).
- When the pumping time is over, the pump will stop and the red warning light will go off.
- Then the sample holder should be inserted completely: with the left hand press lightly against the purple surface on the Compustage and with the right hand turn the sample holder carefully counterclockwise so that the inner pin comes to the 12 o'clock position. The large pin visible on the outside, near the black end of the holder, is thus moved to the 6 o'clock position. This "unlocks" the Compustage. ***Caution: The vacuum pulls the specimen holder in! Therefore, hold it firmly!***
- Slowly slide the specimen holder in to the end. The large pin disappears into the shaft at 6 o'clock. Check whether the specimen grip is completely in.
- Knock loosely against the end of the specimen holder a few times.
- Attach the large gray protective cap.

6. Control Panel

Left Control Panel = LCP

Right Control Panel = RPC

7. Alignments

- **Darken**


Usually the light is off when microscoping. With the dimmer (to the right of the TEM) you can make yourself some dim light.

- **Open Column Valves**

The phosphorescent screen in the TEM must be closed. When the pressure in **Column IGP1** is ≤ 15 , the Column Valves may be opened. The pressure should continue to drop to 6. If the pressure does not reach 20 or 6, then the sample holder should be ejected again and the O-ring inspected for dirt.


- **Load FEG-Registers (Optional)**

If the beam alignment looks very bad, then an older, known FEG register for TEM mode should be loaded. Or of course, if you want to change the mode (TEM, STEM, tomography). To do this, open Setup -> FEG Registers, see figure. Select a register and press Set to load it. The TEM will continue to work with the individual settings of each user where they finished the last session. So if you continue to work in the same mode, then loading new FEG registers should be the exception.

- **Control Spotsize**

At this point you should check whether **Spotsize = 3** (for TEM mode) is set (otherwise change it). Without specimen holder, the TEM automatically switches to Spotsize = 5 for safety reasons.


- **Find the beam**

- Often the beam is on the grid of the sample -> move joysticks a bit (change xy-position). This almost always works. Otherwise you can still try:
- The magnification is too large -> set to 4400 x to 1500 x
- The beam brightness is very far away from the focus point -> change the intensity at the LCP
- Maybe the beam is not directed to the center of the screen -> make magnification very small and change the position with the trackball at the LPC
- If nothing works: Load another FEG register.

- **Condenser Aperture C2 centered**

- For TEM mode, C2 is set to 3.
- Increase magnification slightly, at least 44,000 x.
- Use Intensity on the LCP to focus the beam to a small spot.
- Use the trackball on the LCP to center the spot beam absolutely on the phosphorescent screen.

Faster method

- Then wind the beam fully clockwise with Intensity so that almost the entire screen is illuminated.
- Use the X and Y adjustment knobs of the aperture to adjust it so that the expanded beam is as centered as possible. The four outer markings on the phosphorescent screen can be used as a guide.
- The whole procedure is iteratively repeated until the beam remains centered and uniformly wound (both clockwise and counterclockwise).

More accurate method

- Then pull up the beam with Intensity clockwise until the 4 cm ring on the screen is illuminated.
- Place the trackball exactly on the 4 cm ring.
- Now move Intensity into subfocus (counterclockwise) and illuminate the 4 cm ring again. The difference between the 4 cm ring and the round beam is the misalignment.
- Now using the X and Y adjustment knobs, adjust the aperture so that the beam is slightly better to the 4 cm ring (but not perfect).
- And back to overfocus with Intensity and see how well it fits now.
- Repeat the whole procedure iteratively until the beam stays centered and winds up uniformly (both clockwise and counterclockwise).

• Set eucentric height

There are several methods, all of which should be mastered. Preferably, the sample is always shifted in the z-direction so that it is at the eucentric height (rather than shifting the focus to the sample). This, among other things, makes the calibration of the scale display most accurate, since all base alignments of the TEM were made for the eucentric height.

Method: Minimum contrast

- Press Eucentric Focus on the RCP. This sets the focus for the eucentric height defined in the base alignments of the instrument.
- Set a magnification of 30,000 - 50,000 x.
- Searching for a specimen detail
- Use the z-axis buttons on the RCP to move the specimen up and down so that the specimen detail has minimal contrast.

Method: Minimum scattering


- Press Eucentric Focus on the RCP.
- Set a magnification of 30.000 - 50.000 x.
- Find a sample detail. Collections of amorphous material are best.
- Condense the beam to a point using Intensity on the LCP.
- Use the z-axis buttons on the RCP to move the specimen up and down so that the scattering around the central bright spot is minimized.

Method: focus wobbler

- Press Eucentric Focus on the RCP.
- Set a magnification of 30,000 - 50,000 x. (Also depends strongly on the size of the specimen detail).
- Search for a specimen detail.
- Now fix the sample detail well with the eyes and start the focus wobbler with the "Wobbler" button on the RCP. This will cause the beam to wobble around the specimen (similar to the α -wobbler where the specimen wobbles in the fixed beam - the relative motion is the same).
- Amplitude and spatial direction of the wobble can be changed with the two multifunction controllers.
- Wobbling produces two superimposed images of the sample. Use the z-axis buttons on the RCP to move the specimen up and down so that the two images overlap again. If necessary, correct with the joystick so that the specimen detail does not move out of the image.
- If the eucentric height is correct, then you can set the wobbler amplitude both high and low without changing the image. Turn the wobbler off again.

Method: α -Wobbler

- Press Eucentric Focus on the RCP.
- Set a magnification of 30,000 - 50,000 x. (This also depends on the size of the specimen detail).
- Find a sample detail and move it to the center of the screen.
- Now start the α -wobbler. This causes the specimen holder to tilt periodically about the α -axis. The α -wobbler can be controlled either via one of the multifunction keys L1 or via the user interface of the Compustage (see figure). There you can also set the amplitude, e.g. to 5.



- The z-axis buttons on the RCP are used to move the specimen up and down so that the movement of the specimen detail is minimized.
- When the movement of the specimen detail is minimal, then one is in focus and the α -wobbler can be issued.

- **Alignment: Gun Tilt (Optional)**

Gun Tilt and Gun Shift only need to be redone if the beam intensity has become too weak. **Typical values for the Screen Current ~10 nA and Exposure Time ~0.5 s.**


- Set typical magnification, e.g. 44,000 x
- Illuminate the 4 cm ring (wound up clockwise)
- Click Alignments -> Direct Alignments -> Gun Tilt.
- Now use the two multifunction controls to move the beam back and forth so that the Exposure Time is set to minimum. If necessary, use the trackball to bring the beam back to the center.
- End with Done.

- **Alignment: Gun Shift (Optional)**


Gun Tilt and Gun Shift only need to be redone if the beam intensity has become too weak.

- Set a typical magnification, e.g. 44.000 x.
If the spot moves out of the field of view, start with a smaller magnification.
- Set Spot Size to 9 with R3.
- Use Intensity to contract the beam and center it with the trackball.
- Set Spotsizes to 3 with L3. (Now it is no longer in the center if Gun Shift is not ideal).
- Click Gun Shift in the Direct Alignments.
- Spot the beam with Intensity and center it with the two multifunction sliders.
- Press Done when done.
- Repeat the last steps iteratively until the beam stops moving when the Spotsizes are changed from 9 to 3.
- When Exposure Time has increased again, then readjust Gun Tilt again.
- If necessary, correct the astigmatism of the Condenser Lens again.



- **Alignment: Astigmatismus der Condenser Lens**

- If the beam on the phosphorescent screen is not round, the astigmatism of the Condenser Lens must be corrected.
- Press the stigmator on the LCP.
- The correct stigmator (here: Condenser) must be selected on the PC interface. The menu is at Alignments -> Stigmator.
- With the two multifunctional controllers, the beam can now easily be brought into a round shape. The easiest way is to orientate at the 4 cm ring on the phosphorescent screen.
- Finally, press Stigmator again to turn off the control.

- **Alignment: Pivot Points**

- Set typical magnification, e.g. 175,000 x
- Use Intensity to constrict the beam to a point and center it with the trackball.
- Click on Beam tilt pp X. The beam now wobbles back and forth.
- Adjust the beam with the two multifunction controls so that it has a bright point in the center around which the beam wobbles (see illustration).
- Press Done when done.
- Do the same for Beam tilt pp Y.

Pivot Points: Correct


Pivot Points: Wrong

Pivot Points: Wrong


- **Alignment: Beam Shift**

- Set the typical magnification, e.g. 175,000 x.
- Use Intensity to contract the beam to a point and center it with the trackball.
- Press Beam Shift. The beam now moves to the side. If it disappears from the field of view, reduce the magnification until you find it again.
- Center the beam with the two multifunction controls.
- Press Done when finished.

- **Alignment: Rotation Center**


- Set the typical magnification, e.g. 175,000 x.
- Use Intensity to contract the beam to a point and center it with the trackball.
- Press Rotation Center. The beam now begins to oscillate. Use the ring on the focus control to change the amplitude of the oscillation.
- With the two multifunction knobs you can influence the movement of the beam. It should be adjusted so that it oscillates in the center (and not back and forth), see figure.
- Press Done when finished.
- If Rotation Center is set, then briefly check Beam Shift and Rotation center again, since the alignments influence each other.

Rotation Center: Correct
The beam oscillates
around a bright spot in
the center

Rotation Center: Wrong
The beam oscillates
around a bright point at
the edge

Rotation Center: Quite wrong
The beam oscillates back and forth
so much that a dumbbell-shaped im-
age is created

- **Alignment: Coma-Free (Optional)**


Only necessary at maximum resolution for best images. Coma-Free Alignment and Rotation Center influence exactly the same beam parameters - only the adjustment method is different.

- Set Magnification very large, 440,000 x
- Place a sample detail in the center of the screen and pull the beam wide open with the Intensity.
- Now press Coma-Free Alignment X. This will cause the beam to start bouncing back and forth. The beam should be pulled up so far that the specimen detail always remains visible.
- Now use the multifunction knob on the LPC to focus the specimen detail as it jumps back and forth. The best way to see this is in Digital Micrograph in the FFT. The sample detail itself will look slightly different as it bounces back and forth, because the beam is shining through it in different ways.
- The same for Coma-Free Alignment Y.


- Alignment: Astigmatism of the Objective Lens**

This alignment is made with the support of FFT in Digital Micrograph.

- Holding down the ALT key will bring up a dotted square on the live image. (Start far left at the top and then drag it up).
- Then the Fast Fourier Transform is taken from it.
- Process -> Live -> Reduced FFT

- In the live image of the FFT you can see circles if the eucentric height is not set quite correctly.
- Now you can adjust the eucentric height ideally with z-axis at the RCP. Then the image is well illuminated in the FFT.
- Then use the focus knob to go back to subfocus a bit to see the circles in the FFT.


- If the circles in the FFT are oval, then the astigmatism of the objective lens must be corrected. To do this, press the stigmator on the LCP. If the astigmatism is very strong, the objective lens aperture is probably not set well (correct again if necessary).
- The correct stigmator (here: Objective) must be selected on the PC interface. The menu is at Alignments -> Stigmator. By default, the Objective stigmator is already selected when the Stigmator button is pressed.
- Now use the two multifunction controls to make the circles as round as possible. Sometimes it is helpful to place a circle from the upper menu in the FFT image to see exactly when the circles are round.
- Finally, press Stigmator again to make out the control.
- Use the focus button to go back into focus.


8. Measuring with Digital Micrograph

- **Save images:**


- With "Save Numbered" you only have to enter the sample designation and the date once. Each time you save, the images are automatically labeled and assigned a sequential number.

- File -> Global -> Global Info -> Saving -> Save Numbered
- Then specify there:
 - The path where it should be saved (please, the own directory).
 - The string that should identify the files. Typically sample name + date.
 - Set Next Index to 1.
 - Save as Gatan Format.
 - Confirm with OK.


- **Capture images:**
 - The beam should illuminate the 4 cm circle well.


- Open the phosphorescent screen so that the electron beam falls on the CCD camera. To do this, press either R1 on the RCP or the mechanical button directly on the left of the TEM.

Mechanical button
for lift screen


- With the "Start/Stop View" button in the right menu the camera takes live pictures in fast search mode. (If the camera is still retracted at the beginning, DM asks if you want to extend the camera: Yes).
- Camera modes:
 - Search: Lower resolution (with binning 2) with fast repetition rate for real-time images.
 - Preview: Higher resolution with medium repetition rate.
 - Acquire: Acquire the image with the set exposure time.
- A dotted square is placed on the live image. To do this, select the corresponding symbol in the menu, hold down the ALT key, click in the upper left corner of the image and then drag to the lower right. Then the Fast Fourier Transform is created from it. Under Process ->Live -> Reduced FFT. The FFT is used to control focus and astigmatism (see Alignments).
- Normally you only need Search and Acquire mode. Both are set to the same 2004 x 1336 field, where Search works with 0.1 seconds and binning 2, Acquire with binning 1. You can check this by pressing the small tool icon.
- Normally you only set the exposure time e.g. 1 second.


9. Contrast Aperture

A better contrast is obtained with the objective aperture as the contrast aperture. However, an aperture also limits the maximum resolving power that can be achieved. Typical values:
At 50,000 x aperture 2, at high resolution 350,000 x aperture 4 or no aperture.

- The phosphorescent screen must be closed. Otherwise the scintillator of the camera may be damaged! The objective aperture is still out.
- The sample site that is currently visible in the TEM will be strongly irradiated in the next step. It is therefore advisable to move the joystick slightly to the side to avoid damaging the area of interest.
- Diffraction button on the RCP. This now maps the intermediate image plane of the Objective Lens onto the phosphorescent screen. The spot in the center is very bright. The "Magnification" button is used to set the camera length in Diffraction Mode; set 460 mm here. The spot in the center can be adjusted with the two multifunction controls.



- Now use the control dial to set the aperture size and bring in the aperture (lever to the left).
- Using the two x and y adjustment knobs on the aperture, the aperture must now be set so that the direct beam passes through as centrally as possible.
- For aperture 1-2, the apertures are very narrow compared to the beam. Therefore, it is best to move the small auxiliary screen in with the lever on the left of the TEM and look at the beam through the optical microscope when adjusting the aperture.
- Finally, press the Diffraction button again and return to Bright Field.


10. Ending the Session

- Set magnification to SA 4400 x and illuminate 4 cm circle on screen
- Drive Compustage to the center
First move the joystick to the center, then the rest can be done with Compusage -> Stage² Flapout -> Reset: Holder.
- Take out objective & SA apertures
I. e. both levers to the right.
- Close Column Valves

- Stop camera and move back
"Stop View" and uncheck "Camera Inserted" in Digital Micrograph.
- Remove the sample holder from the Compustage
 - Remove the cover of the Compustage.
 - Touch the Selected Area aperture control once with your left hand to discharge any static charges. Now stand still.
 - Then press against the purple surface on the Compustage with your left hand. This is important, otherwise the Compustage can be damaged when the sample holder is pulled out.
 - **Carefully** pull out the sample holder until you feel a slight resistance. Then turn the sample holder clockwise by about 180° until it stops moving. This "closes" the compustage.
 - For the double tilt holder: Disconnect cable connection.
 - Now carefully pull out the sample holder against the resistance. Stabilise with your left hand. At a certain point it is very easy. Then do not hit the sample holder against the opening.
 - If the red warning light is still illuminated on the Compustage: Select "no holder" on the control interface on the PC and confirm with the small arrow button.
 - Finally, fit the cover of the Compustage.
- Remove sample
 - The part between the O-ring and the head end of the holder must never be touched by hand or other non-vacuum-compatible parts. Wear protective gloves if necessary.
 - Return the sample holder to the holding device.
 - Insert the tool into the hole at the head end of the clamp. Then **carefully** lift the clamp straight up **without force** until it is at right angles to the holder.
 - Remove the sample. Trick: Turn the sample holder and drop the sample directly into the sample box.
 - Close the clamp again.


- **Refit the empty sample holder in the Compustage**
 - This keeps the holder clean and reduces the amount of dirt that gets into the Compustage
 - Finally, slide the protective cap onto the holder.

- **Set FEG to 4500 V (if nobody wants to measure afterwards)**

- **Pour back the remains of liquid nitrogen (if nobody wants to measure afterwards)**
 - Cover the viewing screen on the TEM with the rubber cover.
 - Put on safety goggles and protective gloves.
 - Pour the remainder of the cold trap and the transfer container back into the large nitrogen container.
 - The cold trap and the transfer container are placed upside down in the polystyrene device so that condensation can drop out of the containers.
 - Place a paper towel where the cold trap was on the TEM to collect condensation water.

- **Start Cryo Cycle (if nobody wants to measure afterwards)**
 1. Setup –Flapout > Cryo ->Press Cryo Cycle.
 2. the Cryo Cycle starts automatically after 7 minutes. The TEM user interface can be closed, the programme for the Cryo Cycle continues to run in the background.


- **Exit the three programmes (File -> Exit)**

- **Log out of the TEM control computer (DO NOT shut down), switch off the monitors**
 - **Attention:** A short waiting time (at least 10 seconds) should be observed between closing the three programmes and logging out so that all programmes are really closed.
 - **Attention:** Sometimes the TEM control computer freezes or crashes immediately after logging out. Then inform the TEM instrument officer.

- **Move data**
Data from TEM Bilder:\Nutzer\Username move to TEM:\Nutzer Verzeichnisse\Username.

- **Batch Convert of the DM3-Data**

The images are in DM3 format and each 42 MB in size. They can be edited with Digital Micrograph or with the freeware imageJ. To convert the images in Digital Micrograph into common formats: File -> Batch Convert. Then select the folder where the images are saved and select the format to which you want to convert (BMP, JPEG or TIFF). Confirm with OK and the converted images will be created in the same folder.

- **Make any final entries in the logbook and shut down the computer in the anteroom (if nobody wants to measure afterwards)**

- **Clean up!**


The laboratory should be left at least as tidy and clean as you found it.

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
12. Selected Area Diffraction


The diffraction mode images the back focal plane onto the screen or detector. The diffraction pattern is an image of the electrodes diffracted and scattered in the sample. ATTENTION: The single points (especially the point in the center) of the diffraction pattern can become very intense and thus damage the scintillator of the CCD camera! Therefore:

- **This imaging mode is only allowed for specially trained users!**
- **always work with the utmost caution!**

- Mount sample and insert (for crystalline samples take double tilt holder if necessary, if α AND β are to be tilted).
- For double tilt holder: Sample, washer, conical ring. First turn a little counterclockwise so that you find the place where the thread goes in, then turn clockwise very carefully without force a little pure and "butter soft" only tighten. It must be flush (look at it from the side). For cross section specimens: Install the adhesive edge perpendicular to the rod.
- Mount specimen as in normal TEM mode (bright field), set C1, set eucentric height, perform alignments directly.
- Align specimen with α and β : Press Diffraction (screen must be closed, of course!). In the image, if the sample is thick enough, you can also see Kikuchi lines. Then tilt with α and β , so that in the pole. Finally press Diffraction to go back to bright field.

DP with Kikuchi lines.


By tilting, the sample is now in the pole.

The photos show the 111 direction of silicon.

Examples of Kikuchi lines in various materials can be found in Wiliams Carter Volume 2 (available in the TEM anteroom). Diffraction patterns and Kikuchi lines for various materials can be calculated online at: <http://emaps.mrl.uiuc.edu/>

- Bring the interesting specimen detail to the center. Set the selected area aperture (lowest aperture), set the aperture size (usually no. 2 or even 1) and center it with the two sliders.
- Screen closed to protect camera! Now press the Diffraction button on the RPC. Now the rear focal plane is imaged.
- Use the multifunction knobs to center the centrally transmitted beam. Use the magnification knob to set the camera length, which corresponds to magnification in this figure (typically: 285 or 460 mm).
- Spread the beam with the "Intensity" control on the LPC until the pattern shows rings, ideally with dots for crystalline samples.
- Focus the diffraction pattern with the focus button RCP. If necessary, you can also see this well with the light microscope on the auxiliary screen

- Bring in beam stopper so that the central spot is blocked. This is very important to protect the camera!

On the right of the TEM column is the control for the beam stopper

Turn here to bring the beam stopper in/out

Turn here to change the length of the beam stopper slightly.

- Possibly readjust the central spot with the multifunction controllers so that it is blocked by the beam stopper.



- In the following: DO NOT use the search mode in Digital Micrograph (or TIA), otherwise the camera will be illuminated during each search process. If necessary, press "Stop-View" in DM. To be on the safe side, always close the screen quickly after each image.

- In Digital Micrograph: Acquire with Exposure time 0.01 s (with Binning 1, which is set by default). Acquire an image, close the screen and read the max. intensity of this image. To do this, click on the captured image and call up the statistics under Analysis->Statistics->all and read off "Frame: Maximum...".


The max intensity should not exceed 7.000. Now take some test images with higher times (up to 1 s). If the intensity and the expo-sure time are well adjusted (good image at not more than 7,000 max. intensity), then take and save the image.

-

- Ready. Screen closed. Exit Diffraction mode. Beam stopper out. SA aperture out.

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